NADAL[®] COVID-19 Ag+Influenza A/B Test (test cassette)

REF 243203N-20

DE Gebrauchsanweisung	2	CZ Návod k použití	44
EN Instructions for use	8	FI Käyttöohje	50
FR Instructions d'utilisation	14	SE Användarinstruktioner	56
ES Instrucciones de uso	20	DK Brugervejledning	62
IT Istruzioni per l'uso	26	NO Bruksanvisning	68
PL Sposób użycia	32	NL Gebruiksaanwijzing	74
PT Instruções de Utilização	38	Symbols	83
		Our Teams	84

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Rev. 0, 2020-11-10 OM

1. Intended Use

The NADAL® COVID-19 Ag+Influenza A/B Test is a lateral flow chromatographic immunoassay for the simultaneous qualitative detection of SARS-CoV-2 as well as influenza virus type A and B viral nucleoprotein antigens in human nasopharyngeal specimens (see section 12 'Limitations'). This test is intended for use as an aid in the diagnosis of infections with SARS-CoV-2 as well as influenza virus type A and B. Note that the concentration of viral nucleoprotein antigens may vary in the course of the disease and might fall below the detection limit of the test. Possible infectiousness of test subjects cannot be ruled out based on negative test results. The test procedure is not automated and requires no special training or qualification. The NADAL® COVID-19 Ag+Influenza A/B Test is designed for professional use only.

2. Introduction and Clinical Significance

COVID-19 (Corona Virus Disease) is the infectious disease caused by the recently discovered coronavirus SARS-CoV-2. The most common symptoms of COVID-19 are fever, dry cough, fatigue, sputum production, shortness of breath, sore throat and headache. Some patients may have myalgia, chills, nausea, nasal congestion and diarrhoea. These symptoms begin gradually and are mild in most of the cases. Some people become infected but do not develop any symptoms and do not feel unwell. Most people (about 80%) recover from the disease without special treatment. Approximately one in six people who get infected with COVID-19 becomes seriously ill and develops difficulty breathing. Elderly people, and those with pre-existing conditions, such as high blood pressure, heart problems or diabetes, are more likely to develop serious illness. COVID-19 is transmitted via respiratory droplets that are exhaled by infected people via coughing, sneezing or talking. Most estimates of the incubation period for COVID-19 range from 1 to 14 days, during which people might already be infectious without showing disease symptoms.

Influenza is a highly contagious viral infection of the upper respiratory tract which is characterised by antigen variability, seasonality and the impact on the general population. Of the two main types (A and B) of influenza viruses, influenza A subtypes are differentiated by the antigen variability of surface glycoproteins (haemagglutinin and neuraminidase). The influenza A virus is the most prevalent and is associated with the most serious epidemics. Influenza can cause severe complications such as bronchitis or pneumonia, particularly in children, the elderly or those with chronic respiratory diseases. However, it most commonly occurs as a mild viral infection transmitted by respiratory secretions through sneezing or coughing. There are many other viral infections that can mimic influenza symptoms, making laboratory tests necessary to distinguish it from other acute respiratory infections.

With a sensitivity and specificity of almost 100%, RT-PCR is the gold standard of laboratory diagnostics. However, as a laboratory method, it does not allow immediate diagnosis and requires qualified personnel and expensive laboratory equipment. Antiviral influenza drugs, which have been available since the late 1990s, are most effective if administered soon after the onset of the illness (within 48 hours). Rapid tests for the detection of influenza antigens can

therefore enable successful antiviral therapy to be carried out in good time, thereby improving patient healthcare and reducing health costs.

3. Test Principle

The NADAL® COVID-19 Ag+Influenza A/B Test is a lateral flow chromatographic immunoassay for the qualitative detection of SARS-CoV-2 as well as influenza virus type A and B viral nucleoprotein antigens in human nasopharyngeal specimens.

The NADAL® COVID-19 Ag+Influenza A/B Test contains two internal test strips (left side of the test cassette: COVID-19 Ag Test; right side of the test cassette: Influenza A/B Test).

Anti-SARS-CoV-2 antibodies are immobilised in the test line region 'T' of the membrane. Anti-influenza virus type A and B antibodies are immobilised in the test line regions 'A' and 'B', respectively, of the membrane. A specimen is added to an extraction tube containing buffer in order to release SARS-CoV-2 as well as influenza virus type A and B viral nucleoprotein antigens. During the test, extracted antigens, if present in the specimen, bind to anti-SARS-CoV-2 antibodies and/or anti-influenza virus type A and/or B antibodies conjugated to coloured particles and precoated onto the sample pad of the test cassette. The mixture then migrates along the membrane chromatographically by capillary action and interacts with the reagents on the membrane. The complexes are then captured by anti-SARS-CoV-2 antibodies as well as by anti-influenza virus type A and B antibodies in the test line regions 'T' as well as 'A' and 'B', respectively. Excess coloured particles are captured in each control line region 'C'. The presence of a coloured line in the test line region 'T' and/or 'A' and/or 'B' indicates a positive result. The absence of a coloured line in the test line region 'T' and/or 'A' and/or 'B' indicates a negative result.

The formation of a coloured line in each control line region 'C' serves as a procedural control, indicating that the proper volume of specimen has been added and membrane wicking has occurred.

4. Reagents and Materials Supplied

- 20 NADAL® COVID-19 Ag+Influenza A/B test cassettes*
- Additional material provided according to 93/42/EEC:
 - Due to possible supply shortages of COVID-19 related accessory medical products, the swab manufacturer might change. Therefore, the supplied swabs are from one of the manufacturers listed below.
 - a) 20 sterile swabs, CE0086
 -  Puritan Medical Products Company LLC
31 School Street
Guilford, Maine 04443-0149 USA (authorised EU representative EMERGO EUROPE, The Hague, The Netherlands)
 - b) 20 sterile swabs, CE0197
 -  Jiangsu Changfeng Medical Industry Co., Ltd
Touqiao Town, Guangling District, Yangzhou,
Jiangsu 225109 China (authorised EU representative Llins Service & Consulting GmbH, Obere Seegasse 34/2, 69124 Heidelberg, Germany)

c) 20 sterile swabs, CE0197

CITOTEST LABWARE MANUFACTURING CO., LTD
No.48, Xinxiu Road, Haimen, Jiangsu province
(authorised EU representative WellKang Ltd,
16 Castle St, Dover, CT16 1PW, UK)

d) 20 sterile swabs, Copan Floqswabs; CE0123

Copan Italia S.p.A.,
Via Perotti 10, 25125 Brescia, Italy

e) 20 sterile swabs, CE0123

Zhejiang Gongdong Medical Technology Co., Ltd
No.10 Beiyuan Ave., Huangyan, 318020 Taizhou,
Zhejiang, People's Republic of China (authorised
EU representative Shanghai International Holding
Corp. GmbH (Europe), Eiffestraße 80, 20537
Hamburg, Germany)

- 20 extraction tubes incl. dropper caps
- Buffer (Due to possible supply shortages of COVID-19 related accessory medical products, buffer is supplied either in 2 bottles (7 mL each) or 1 bottle (10 mL). Both options are sufficient for 20 determinations.)**
- 1 reagent holder
- 1 package insert

*containing the preservative sodium azide: <0.1%

**Buffer containing the following preservatives: sodium azide: <0.1 mg/mL
The detergents contained in the buffer lyse and neutralise the virus.

No hazard labelling is required according to Regulation (EC) № 1272/2008 CLP. Concentrations are below exemption threshold.

5. Additional Materials Required

- Timer

6. Storage & Stability

Test kits should be stored at 2-30°C until the indicated expiry date. Test cassettes are stable until the expiry date printed on the foil pouches. Test cassettes must remain in the sealed foil pouches until use. Do not freeze the test kit. Do not use tests beyond the expiry date indicated on the packaging. Care should be taken to protect test kit components from contamination. Do not use test kit components if there is evidence of microbial contamination or precipitation. Biological contamination of dispensing equipment, containers or reagents can lead to inaccurate results.

7. Warnings and Precautions

- For professional *in-vitro* diagnostic use only.
- Carefully read through the test procedure prior to testing.
- Do not use the test beyond the expiration date indicated on the packaging.
- Do not use test kit components if the primary packaging is damaged.
- Tests are for single use only.
- Do not add specimens to the reaction area (result area).
- In order to avoid contamination, do not touch the reaction area (result area).

- Avoid cross-contamination of specimens by using a new extraction tube for each specimen obtained.
- Do not substitute or mix components from different test kits.
- Do not use the buffer if it is discoloured or turbid. Discolouration or turbidity may be a sign of microbial contamination.
- Do not eat, drink or smoke in the area where specimens and test kits are handled.
- Wear appropriate personal protective equipment, such as face mask, isolation gown, gloves and eye protection during specimen collection, preparation and test procedure.
- Handle all specimens as if they contain infectious agents. Observe established precautions for microbiological risks throughout all procedures and standard guidelines for the appropriate disposal of specimens.
- Further specimen processing and patient management should follow local guidelines and regulations.
- The test kit contains products of animal origin. Certified knowledge of the origin and/or sanitary state of the animals does not completely guarantee the absence of transmissible pathogenic agents. It is therefore recommended that these products be treated as potentially infectious and handled in accordance with usual safety precautions (e.g., do not ingest or inhale).
- Temperature can adversely affect test results.
- Used testing materials should be disposed of according to local regulations.

8. Specimen Collection and Preparation

Nasopharyngeal specimen:

- Insert the swab into the nostril, parallel to the palate (not upwards) until resistance is encountered or the distance is equivalent to that from the ear to the nostril of the patient, indicating contact with the nasopharynx.
- Gently rub and roll the swab. Leave swab in place for several seconds to absorb secretions.
- Slowly remove the swab while rotating it. Specimens can be collected from both nostrils using the same swab, but it is not necessary to collect specimens from both sides if the tip is saturated with fluid from the first collection.

Note:


- Use only synthetic fibre swabs with plastic shafts. Do not use calcium alginate swabs or swabs with wooden shafts, as they may contain substances that inactivate some viruses and inhibit further testing.
- Swab specimens should be tested immediately after collection. Use freshly collected specimens for best test performance.
- If not tested immediately, swab specimens can be stored in viral transport media (VTM) without denaturing agents at 2-8°C for 24 hours after collection.
- Viral transport media (VTM) without denaturing agents as well as VTM containing virus-inactivating agents and evaluated with the test can be used during specimen storage for the subsequent antigen detection using the NADAL® COVID-19 Ag+Influenza A/B Test. In order to influence the sensitivity as little as possible, a low volume of VTM (max. 1 mL). The test procedure when using viral transport media (VTM) is available separately:

- under the following link:
<https://www.nal-vonminden.com/en/vtm-ifu>
 - by request at +49 941 29010-0 or
info@nal-vonminden.com
- Do not use specimens that are obviously contaminated with blood, as it may interfere with the flow of specimens and lead to inaccurate test results.



9. Test Procedure

Bring tests, specimens, buffer and/or controls to room temperature (15-30°C) prior to testing.


1. Place a clean extraction tube labeled with the patient or control identification into the designated area of the reagent holder.
2. Gently mix the buffer by carefully swiveling the bottle.
3. Holding the buffer bottle vertically and without touching the edge of the tube, add 10 drops to the extraction tube.
4. Insert the swab with the collected specimen into the tube. Swirl the swab and squeeze it 10-15 times by compressing the wall of the extraction tube against the swab to extract the antigens contained in the swab.

5. Let the solution stand for 2 minutes.
6. Remove the swab, pressing it firmly against the wall of the tube to release as much liquid as possible. Dispose of the swab in accordance with guidelines for the handling of infectious agents.
7. Remove the test cassette from the foil pouch and use it as soon as possible. The best results will be obtained if the test is performed immediately after opening the foil pouch. Label the test cassette with the patient or control identification.

8. Place the test cassette on a clean and level surface.
9. Attach a dropper cap to the extraction tube, invert the tube and transfer 2 drops of the extracted solution to each specimen well (S) of the test cassette.

10. Start the timer.

11. Wait for the coloured line(s) to appear. Read the test result after exactly 15 minutes. Do not interpret the result after more than 15 minutes.

10. Result Interpretation


Positive for COVID-19 Ag:

A coloured line develops in the control line region 'C' and another coloured line develops in the test line region 'T'.

Positive for influenza virus type A:

A coloured line develops in the control line region 'C' and another coloured line develops in the test line region 'A'.

Positive for influenza virus type B:

A coloured line develops in the control line region 'C' and another coloured line develops in the test line region 'B'.


Note: A combination of positive results for various parameters is possible.

Note: The colour intensity in the test line region 'T', 'A' or 'B' may vary depending on the concentration of SARS-CoV-2 as well as influenza virus type A and B viral nucleoprotein antigens in the specimen. Any shade of colour in the test line region 'T', 'A' or 'B' should be considered a positive result. Note that this is a qualitative test only and it cannot determine the analyte concentration in the specimen.

Negative:

A coloured line develops in each control line region 'C'. No lines develop in the test line regions 'T', 'A' and 'B'.

Invalid

The control line 'C' fails to appear. Results from any test which has not produced a control line at the specified reading time must be discarded. Please review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your distributor.

Insufficient specimen volume, incorrect operating procedure or expired tests are the most likely reasons for control line failure.

11. Quality Control

An internal procedural control is included in the test cassette: A coloured line developing in each control line region 'C' is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

Good laboratory practice (GLP) recommends the use of external control materials to ensure proper test kit performance.

12. Limitations

- The NADAL® COVID-19 Ag+Influenza A/B Test is for professional *in-vitro* diagnostic use only. The test should be used for the qualitative detection of SARS-CoV-2 as well as influenza virus type A and B viral nucleoprotein antigens in human nasopharyngeal specimens only. Neither the quantitative value nor the rate of increase/decrease in the concentration of SARS-CoV-2 as well as influenza virus type

A and B viral nucleoprotein antigens can be determined with this qualitative test.

- The NADAL® COVID-19 Ag+Influenza A/B Test only detects the presence of SARS-CoV-2 as well as influenza virus type A and B viral nucleoprotein antigens in specimens and should not be used as the sole criterion for a diagnosis of COVID-19 and influenza A or B.
- Both viable and non-viable SARS-CoV-2 viruses as well as influenza viruses type A and B can be detected using the NADAL® COVID-19 Ag+Influenza A/B Test.
- The sections 'Specimen Collection and Preparation' as well as 'Test Procedure' must be followed closely while testing. Failure to follow them may lead to inaccurate test results because the antigen concentration in the swab is highly dependent on the correct procedure.
- As with all diagnostic tests, all results should be interpreted in conjunction with other clinical information available to the physician.
- In the course of SARS-CoV-2 or influenza A or B infection, the concentration of viral nucleoprotein antigens may fall below the detection limit of the test.
- If the test result is negative and clinical symptoms persist, additional testing using other clinical methods is recommended. A negative result does not at any time preclude the possibility of a SARS-CoV-2 or influenza A or B infection and should be confirmed via molecular assay.
- Additional testing is required to differentiate any specific influenza virus A and B subtypes or strains in consultation with state or local public health departments.
- Monoclonal antibodies may fail to detect or detect with less sensitivity influenza A virus antigens that have undergone minor amino acid changes in the target epitope region.
- Individuals who received nasally administered influenza A vaccine may have positive test results for up to three days after vaccination.
- Positive results do not preclude co-infections with other pathogens.
- Children tend to shed influenza virus more abundantly and for longer periods of time than adults. Therefore, testing specimens from adults may often yield lower sensitivity than testing specimens from children.
- Positive and negative predictive values are highly dependent on prevalence. The local prevalence should be taken into consideration when interpreting diagnostic test results.

13. Expected Values

SARS-CoV-2 viral particles are normally present in the respiratory tracts of COVID-19 patients. A positive test result can indicate an acute infection. Virus concentrations in nasopharyngeal swab specimens may vary in the course of the disease and might fall below the detection limit of rapid tests, even though patients are still showing symptoms. Conversely, the virus might continue to be detectable over long periods of time even in convalescent patients. Possible infectiousness of test subjects cannot be ruled out based on negative test results. On average, 3 to 10% of unvaccinated adults get infected with influenza each year (the number of infections varies significantly depending on the currently most prevalent strains). A positive rapid test result can indicate the presence of influenza virus type A and/or B viral nucleoprotein antigens in

the specimen material, suggesting an acute infection with influenza A and/or B. Since patients might also acquire co-infections with other severe respiratory diseases, additional testing for other viral and bacterial infectious agents is recommended.

14. Performance Characteristics

Clinical performance (specimens extracted in buffer without VTM)

Diagnostic sensitivity and specificity

The NADAL® COVID-19 Ag+Influenza A/B Test was evaluated with clinical specimens whose status was confirmed using RT-PCR (C_t range positive: 20-37). The sensitivity for SARS-CoV-2 was calculated for the range from high to medium viral load (C_t 20-30) and from high to very low viral load (C_t 20-37).

The results are presented in the following tables.

NADAL® COVID-19 Ag+Influenza A/B Test (COVID-19 Ag)	RT-PCR, C_t 20-30			
		Positive	Negative	Total
	Positive	120	0	120
	Negative	3	161	164
Total	123	161	284	

Diagnostic sensitivity (C_t 20-30): 97.6% (93.1% - 99.2%)*

Overall agreement (C_t 20-30): 98.9% (96.9% - 99.6%)*

Diagnostic specificity: >99.9% (97.7% - 100%)*

*95% confidence interval

NADAL® COVID-19 Ag+Influenza A/B Test (COVID-19 Ag)	RT-PCR, C_t 20-37			
		Positive	Negative	Total
	Positive	150	0	150
	Negative	37	161	198
Total	187	161	348	

Diagnostic sensitivity (C_t 20-37): 80.2% (73.9% - 85.3%)*

Overall agreement (C_t 20-37): 89.4% (85.7% - 92.2%)*

Diagnostic specificity: >99.9% (97.7% - 100%)*

*95% confidence interval

In order to show the strong dependence of the diagnostic sensitivity on the viral load, the following table demonstrates the SARS-CoV-2 sensitivity for different C_t value ranges of the reference PCR:

C_t range	Sensitivity
20 – 25	97.12%
20 – 30	97.56%
20 – 32	96.21%
20 – 35	85.71%
20 – 37	80.21%

Please note that C_t values may vary between different PCR systems at the same virus concentration.

A total of 190 nasopharyngeal swabs from patients with suspected influenza A viral infection were tested using the NADAL® COVID-19 Ag+Influenza A/B Test. Of these specimens, 78 were confirmed to be positive and 112 were confirmed to be negative by cell culture.

The results are presented in the following table:

NADAL® COVID-19 Ag+Influenza A/B Test (Influenza A)		Cell culture		
		Positive	Negative	Total
		Positive	65	9
Negative	13	103	116	
Total	78	112	190	

Positive agreement with cell culture: 83.3% (73.5%-90%)*

Negative agreement with cell culture: 92% (85.4% - 95.7%)*

Total agreement with cell culture: 88.4% (83.1% - 92.2%)*

*95% confidence interval

A total of 190 nasopharyngeal swabs from patients with suspected influenza B viral infection were tested using the NADAL® COVID-19 Ag+Influenza A/B Test. Of these specimens, 85 were confirmed to be positive and 105 were confirmed to be negative by cell culture.

The results are presented in the following table:

NADAL® COVID-19 Ag+Influenza A/B Test (Influenza B)		Cell culture		
		Positive	Negative	Total
		Positive	70	6
Negative	15	99	114	
Total	85	105	190	

Positive agreement with cell culture: 82.4% (72.9% - 89%)*

Negative agreement with cell culture: 94.3% (88.1% - 97.4%)*

Total agreement with cell culture: 88.9% (83.7% - 92.7%)*

*95% confidence interval

Analytical performance

Analytical sensitivity study (detection limit)

The detection limit of the NADAL® COVID-19 Ag+Influenza A/B Test is $2 \times 10^{2.4}$ TCID₅₀/mL and was determined with a SARS-CoV-2 control with a known virus titre.

The detection limit of the NADAL® COVID-19 Ag+Influenza A/B Test is 0.4 ng/mL for recombinant SARS-CoV-2 nucleoprotein.

Concentrations determined as the detection limit for each influenza A/B strain tested are listed below:

Influenza A: A2/Aichi/2/68(H3N2), 2.3×10^3 CEID₅₀/test*

Influenza B: Hong Kong 5/72, 3.5×10^3 CEID₅₀/test*

*CEID₅₀: Chicken Embryo Infectious Dose

Interfering substances

The following substances, normally present in respiratory specimens or artificially introduced into the respiratory tract, were evaluated at the concentrations listed below and showed no interference with the NADAL® COVID-19 Ag+Influenza A/B Test.

Substance	Concentration	Substance	Concentration
3 OTC nasal sprays	10%	Guaiaicol glyceryl ether	20 mg/mL
3 OTC mouthwashes	10%	Mucin	1%
3 OTC sore throat liquids	10%	Mupirocin	250 µg/mL
4-acetamidophenol	10 mg/mL	Oxymetazoline	10 mg/mL

Substance	Concentration	Substance	Concentration
Acetylsalicylic acid	20 mg/mL	Phenylephrine	10 mg/mL
Albuterol	20 mg/mL	Phenylpropanolamine	20 mg/mL
Chlorpheniramine	5 mg/mL	Relenza® (zanamivir)	20 mg/mL
Dexamethasone	5 mg/mL	Rimantadine	500 ng/mL
Dextromethorphan	10 mg/mL	Tamiflu® (oseltamivir)	100 mg/mL
Diphenhydramine	5 mg/mL	Tobramycin	40 mg/mL
Doxylamine succinate	1 mg/mL	Triamcinolone	14 mg/mL
Flunisolide	3 mg/mL		

Cross-reactivity

Specimens spiked with the following pathogens were tested using the NADAL® COVID-19 Ag+Influenza A/B Test (left side of the test cassette: COVID-19 Ag Test):

HCoV-HKU1, HCoV-OC43, HCoV-NL63, HCoV-229E, measles virus, *Streptococcus pneumoniae*, Epstein-Barr virus, *Bordetella parapertussis*, influenza A (H1N1) pdm09, influenza A (H3N2), influenza A (H5N1), influenza A (H7N9), influenza A (H7N7), influenza B Victoria lineage, influenza B Yamagata lineage, *Haemophilus influenzae*, *Candida albicans*, *Mycobacterium tuberculosis*, respiratory syncytial virus, adenovirus, parainfluenza virus type 1, 2, 3, human metapneumovirus, rhinovirus, coxsackievirus type A16, norovirus, mumps virus, *Legionella pneumophila*, *Mycoplasma pneumoniae*, *Chlamydia pneumoniae*, *Streptococcus pyogenes*, *Streptococcus agalactiae*, group C *Streptococcus*, *Staphylococcus aureus*.

No cross-reactivity with the specimens was observed when tested using the NADAL® COVID-19 Ag+Influenza A/B Test (left side of the test cassette: COVID-19 Ag Test).

Specimens spiked with the following pathogens were tested using the NADAL® COVID-19 Ag+Influenza A/B Test (right side of the test cassette: Influenza A/B Test):

Human adenovirus B, C; adenovirus type 10, 18; coxsackievirus type A9, B5; human herpesvirus type 2, 5; herpes simplex virus type 1; human rhinovirus type 2, 14, 16; mumps virus; Murine respirovirus; respiratory syncytial virus; Rubella virus; SARS-CoV-2, human coronavirus OC43; echovirus type 2, 3, 6; measles virus; parainfluenza virus type 2, 3; varicella-zoster virus.

Acinetobacter calcoaceticus, *Bacteroides fragilis*, *Bordetella pertussis*, *Bacillus cereus*, *Bacillus subtilis*, *Bordetella parapertussis*, *Moraxella catarrhalis*, *Chlamydia pneumoniae*, *Corynebacterium diphtheriae*, *Citrobacter freundii*, *Enterobacter cloacae*, *Enterococcus faecalis*, *Escherichia coli*, *Gardnerella vaginalis*, *Haemophilus influenzae*, *Klebsiella oxytoca*, *Klebsiella pneumoniae*, *Lactobacillus casei*, *Lactobacillus plantarum*, *Legionella pneumophila*, *Listeria monocytogenes*, *Mycobacterium avium*, *Mycobacterium intracellulare*, *Mycobacterium tuberculosis*, *Mycoplasma pneumoniae*, *Neisseria meningitidis*, *Neisseria sicca*, *Neisseria subflava*, *Nocardia asteroides*, *Proteus vulgaris*, *Pseudomonas aeruginosa*, *Serratia liquefaciens*, *Staphylococcus aureus*,

Staphylococcus epidermidis, groups A, B, C, F, G *Streptococcus*, *Streptococcus mutans*, *Streptococcus pneumoniae*, *Streptococcus salivarius*, *Streptococcus sanguis*, *Yersinia enterocolitica*.

Precision

Repeatability and reproducibility

Precision for SARS-CoV-2 and Influenza A/B was established by testing 10 replicates of negative, low positive and high positive controls.

Reproducibility for SARS-CoV-2 was established by testing triplicates of negative, low and high positive controls. Testing was performed by 3 operators using 3 independent NADAL® COVID-19 Ag+Influenza A/B test lots at 3 different sites on 5 separate days.

Reproducibility for Influenza A/B was established by testing 10 replicates of negative, low and high positive controls. Testing was performed using 3 independent NADAL® COVID-19 Ag+Influenza A/B test lots.

The NADAL® COVID-19 Ag+Influenza A/B Test demonstrated acceptable repeatability and reproducibility. The negative and positive values were correctly identified >99% of the time.

15. References

1. Cui J, Li F, Shi ZL, Origin and evolution of pathogenic coronaviruses, *Nat Rev Microbiol* 2019; 17:181-192.
2. Su S, Wong G, Shi W, et al, Epidemiology, genetic recombination, and pathogenesis of coronaviruses, *Trends Microbiol* 2016;24:490-502.
3. Weiss SR, Leibowitz JL, Coronavirus pathogenesis, *Adv Virus Res* 2011; 81:85-164.
4. Murphy, B.R., and R.G. Webster, 1996, Orthomyxoviruses, pp.1397-1445. In: Fields, Virology, 3rd edition, B.N. Fields, D.M. Knipe, P.M. Howley, et al. (eds.), Lippincott-Raven, Philadelphia.
5. Biosafety in Microbiological and Biomedical Laboratories, 4th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (1999).
6. CDC website: <http://www.cdc.gov/flu/>
7. Anne Moscona. Neuraminidase Inhibitors for Influenza, 2005. *The New England Journal of Medicine*, 353 (13):1363-1373.

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a) 20 tamponi sterili, CE0086

Puritan Medical Products Company LLC
31 School Street

Guilford, Maine 04443-0149 USA (rappresentante autorizzato EU EMERGO EUROPE, The Hague, The Netherlands)

b) 20 tamponi sterili, CE0197

Jiangsu Changfeng Medical Industry Co., Ltd
Touqiao Town, Guangling District, Yangzhou,

Jiangsu 225109 China (rappresentante autorizzato EU Llins Service & Consulting GmbH, Obere Seegasse 34/2, 69124 Heidelberg, Germany)

c) 20 tamponi sterili, CE0197

CITOTEST LABWARE MANUFACTURING CO., LTD
No.48, Xinxu Road, Haimen, Jiangsu province
(rappresentante autorizzato EU WellKang Ltd, 16 Castle St, Dover, CT16 1PW, UK)

d) 20 tamponi sterili, Copan Floqswabs; CE0123

Copan Italia S.p.A., Via Perotti 10,
25125 Brescia, Italy

e) 20 tamponi sterili, CE0123

Zhejiang Gongdong Medical Technology Co., Ltd
No.10 Beiyuan Ave., Huangyan, 318020 Taizhou,
Zhejiang, People's Republic of China
(rappresentante autorizzato EU Shanghai
International Holding Corp. GmbH (Europe),
Eiffestraße 80, 20537 Hamburg, Germany)

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15. Bibliografia

1. Cui J, Li F, Shi ZL, Origin and evolution of pathogenic coronaviruses, *Nat Rev Microbiol* 2019; 17:181-192.
2. Su S, Wong G, Shi W, et al, Epidemiology, genetic recombination, and pathogenesis of coronaviruses, *Trends Microbiol* 2016;24:490-502.
3. Weiss SR, Leibowitz JL, Coronavirus pathogenesis, *Adv Virus Res* 2011; 81:85-164.

4. Murphy, B.R., and R.G. Webster, 1996, Orthomyxoviruses, pp.1397-1445. In: Fields, Virology, 3rd edition, B.N. Fields, D.M. Knipe, P.M. Howley, et al. (eds.), Lippincott-Raven, Philadelphia.
5. Biosafety in Microbiological and Biomedical Laboratories, 4th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (1999).
6. CDC website: <http://www.cdc.gov/flu/>
7. Anne Moscona. Neuraminidase Inhibitors for Influenza, 2005. *The New England Journal of Medicine*, 353 (13):1363-1373.

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15. Bibliografia

- Cui J, Li F, Shi ZL, Origin and evolution of pathogenic coronaviruses, *Nat Rev Microbiol* 2019; 17:181-192.
- Su S, Wong G, Shi W, et al, Epidemiology, genetic recombination, and pathogenesis of coronaviruses, *Trends Microbiol* 2016;24:490-502.
- Weiss SR, Leibowitz JL, Coronavirus pathogenesis, *Adv Virus Res* 2011; 81:85-164.
- Murphy, B.R., and R.G. Webster, 1996, Orthomyxoviruses, pp.1397-1445. In: Fields, Virology, 3rd edition, B.N. Fields, D.M. Knipe, P.M. Howley, et al. (eds.), Lippincott-Raven, Philadelphia.
- Biosafety in Microbiological and Biomedical Laboratories, 4th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (1999).
- CDC website: <http://www.cdc.gov/flu/>
- Anne Moscona. Neuraminidase Inhibitors for Influenza, 2005. *The New England Journal of Medicine*, 353 (13):1363-1373.

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15. Referências

1. Cui J, Li F, Shi ZL. Origin and evolution of pathogenic coronaviruses, Nat Rev Microbiol 2019; 17:181-192.
2. Su S, Wong G, Shi W, et al. Epidemiology, genetic recombination, and pathogenesis of coronaviruses, Trends Microbiol 2016; 24:490-502.
3. Weiss SR, Leibowitz JL. Coronavirus pathogenesis, Adv Virus Res 2011; 81:85-164.
4. Murphy, B.R., and R.G. Webster, 1996. Orthomyxoviruses, pp.1397-1445. In: Fields, Virology, 3rd edition, B.N. Fields, D.M. Knipe, P.M. Howley, et al. (eds.), Lippincott-Raven, Philadelphia.
5. Biosafety in Microbiological and Biomedical Laboratories, 4th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (1999).
6. CDC website: <http://www.cdc.gov/flu/>
7. Anne Moscona. Neuraminidase inhibitors for Influenza, 2005. The New England Journal of Medicine, 353 (13):1363-1373.

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15. Reference

1. Cui J, Li F, Shi ZL, Origin and evolution of pathogenic coronaviruses, Nat Rev Microbiol 2019; 17:181-192.
2. Su S, Wong G, Shi W, et al, Epidemiology, genetic recombination, and pathogenesis of coronaviruses, TrendsMicrobiol 2016;24:490-502.
3. Weiss SR, Leibowitz JL, Coronavirus pathogenesis, Adv Virus Res 2011; 81:85-164.
4. Murphy, B.R., and R.G. Webster, 1996, Orthomyxoviruses, pp.1397-1445. In: Fields, Virology, 3rd edition, B.N. Fields, D.M. Knipe, P.M. Howley, et al. (eds.), Lippincott-Raven, Philadelphia.
5. Biosafety in Microbiological and Biomedical Laboratories, 4th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (1999).
6. CDC website: <http://www.cdc.gov/flu/>
7. Anne Moscona. Neuraminidase Inhibitors for Influenza, 2005. The New England Journal of Medicine, 353 (13):1363-1373.

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15. Lähteet

- Cui J, Li F, Shi ZL. Origin and evolution of pathogenic coronaviruses, Nat Rev Microbiol 2019; 17:181-192.
- Su S, Wong G, Shi W, et al. Epidemiology, genetic recombination, and pathogenesis of coronaviruses, Trends Microbiol 2016; 24:490-502.
- Weiss SR, Leibowitz JL. Coronavirus pathogenesis, Adv Virus Res 2011; 81:85-164.
- Murphy, B.R., and R.G. Webster, 1996. Orthomyxoviruses, pp.1397-1445. In: Fields, Virology, 3rd edition. B.N. Fields, D.M. Knipe, P.M. Howley, et al. (eds.), Lippincott-Raven, Philadelphia.
- Biosafety in Microbiological and Biomedical Laboratories, 4th Edition. U.S. Department of Health and Human Services, CDC, NIH, Washington, DC (1999).

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